

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 1 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

**Contents**

A.	INTRODUCTION .....	2
B.	EQUIPMENT .....	2
C.	REAGENTS AND SOLUTIONS .....	3
D.	STANDARDS.....	5
E.	SAMPLE PREPARATION.....	6
F.	ANALYTICAL PROCEDURE .....	6
G.	CALCULATIONS .....	8
H.	HAZARD ANALYSIS.....	8
I.	QUALITY ASSURANCE PLAN .....	9
J.	WORKSHEET.....	11
K.	REFERENCES .....	12
L.	APPROVALS AND AUTHORITIES.....	12

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 2 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

**A. INTRODUCTION**

1. Theory

This procedure is used to screen bovine liver and muscle samples for Flunixin residues at concentrations of  $\geq 50$  ppb for liver and  $\geq 10$  ppb for muscle. Blended liver sample and homogenized muscle sample are extracted with a dilute acidic phosphate buffer and analyzed using a commercial Flunixin ELISA kit (Flunixin Enhanced ELISA kit, Neogen Corporation).

2. Applicability

This procedure is applicable to bovine liver at levels  $\geq 50$  ppb and bovine muscle at levels  $\geq 10$  ppb.

**B. EQUIPMENT**

1. Apparatus

Note: Equivalent apparatus may be substituted for those listed below.

- a. Balance, analytical. – Mettler model MT5, sensitive to 0.0001g.
- b. Balance, top loader. – Mettler model PM2000, sensitive to 0.01 g.
- c. Refrigerator, explosion-proof. – Lab-Line Frigid-Cab model 3566.
- d. Freezer, explosion-proof. – Lab-Line Frigid-Cab model 3552.
- e. Vortex mixer. – Lab-Line Super Mixer model 1290.
- f. Shaker (Two Speed), platform. – Eberhard Ultrasonic bath. –VWR Aquasonic model 150T.
- g. Centrifuge (high speed). – Jouan model KR22.
- h. pH meter. – Orion 611, readable to 0.01 pH.
- i. Volumetric dispenser, bottle-top. – Capable of delivering from 1 to 10 mL.
- j. Stirrer/hot plate. – Corning PC-4420.
- k. Micropipettor. – 20 $\mu$ L fixed.
- l. Micropipettors, Variable. – 10-100  $\mu$ L, 20-200  $\mu$ L, and 100-1000  $\mu$ L.
- m. Centrifuge tubes, with caps, for use with high-speed centrifuge. – Nalgene. #3119-0050, 50-mL capacity; Description: Oakridge. Material: polyallomer.
- n. Glass storage bottles with caps. – 60-mL capacity, Daigger # LX2223C and 125-mL capacity, Daigger #LX2223C.
- o. Scintillation vials. – Disposable glass, with tops, Kimble #74505-20.

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 3 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

- p. Racks. – To hold Nalgene centrifuge tubes and working standards stored in scintillation vials.
- q. Volumetric flasks, Class A – 50, 100, 200, 500 and 1000 mL.
- r. Pipettes, Class A. –1, 2, 5, 10 and 25 mL.
- s. Graduated cylinders, Class A. – 100 and 250 mL.
- t. HPLC-grade water system. – Ion pure plus 150.

2. Instrumentation

Note: Equivalent instrumentation may be substituted.

- a. Plate Reader. – Biotec Autoreader, model ELx 808, equipped with 630 nm or 650 nm filter (ELISA Technologies). Reader should be calibrated to assure overall accuracy at the wavelength used.

**C. REAGENTS AND SOLUTIONS**

Note: Equivalent reagents may be substituted.

1. Reagents

- a. Test kits. – Flunixin Enhanced ELISA kit, Cat # 101910 (NEOGEN Corporation, ELISA Technologies Division, 628 E. 3rd Street, Lexington, Kentucky, 40505). Contents of bulk kit include the following:
  - i. EIA buffer (200 mL).
  - ii. 10x Wash buffer concentrate (100 mL). Dilute before use according to manufacturer instructions.
  - iii. K-Blue Substrate (100 mL). Stabilized 3,3',5,5'-Tetramethylbenzidine (TMB) plus hydrogen peroxide in a single bottle. Light sensitive.
  - iv. Drug-Enzyme Conjugate (1 mL). Drug-horseradish peroxidase concentrate. Dilute before use according to manufacturer instructions.
  - v. Pre-coated Plates (5). 96 well Costar plates in strips of 8 breakaway wells coated with anti-drug antisera. Do not wash until the sample/drug-conjugate incubation is complete.

Note: The test kits must be stored in a refrigerator at 2-8°C. Do not use past expiration date.

- b. Sodium phosphate dibasic anhydrous – J.T. Baker # 3828-01, formula  $\text{Na}_2\text{HPO}_4$ , molecular weight 141.96.
- c. Potassium phosphate monobasic crystal – J.T. Baker #-3246-01, formula  $\text{KH}_2\text{PO}_4$ , molecular weight 136.09.

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 4 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

- d. Methanol, HPLC-grade –J.T. Baker # 9093-03.
- e. Water, HPLC-grade – produced from Ion Pure Plus 150 system.
- f. HCl, Concentrated (36.5-38.0%).
- g. Sodium hydroxide pellets – Sigma ACS reagent S-0899. Formula NaOH, Molecular weight 40.00.
- h. Buffer solution, pH 4.0 – Fisher Certified SB98-500.
- i. Buffer solution, pH 7.0 – Fisher Certified SB108-500.

2. Solutions

- a. Extraction Buffer: - Weigh 6.8 g of potassium phosphate monobasic, and 7.1 g of sodium phosphate dibasic into a 1 L class A volumetric flask. Dilute to volume with HPLC grade water. Adjust pH to 6.8 with 1M HCl or 1M NaOH solutions
- b. 1 N HCl solution – Dilute 8.3 mL of concentrated HCl to 100 mL in HPLC-grade water in a 100-mL volumetric flask.
- c. 1 N NaOH solution – Transfer 4 g NaOH into a 100-mL volumetric flask and dilute to mark with HPLC-grade water. Mix thoroughly.
- d. Horseradish peroxidase (HRP) enzyme conjugate - EIA buffer dilution (1+180). (Both reagents are supplied with the bulk kit). Prepare in ratio of 1 µL HRP conjugate to 180 µL of EIA buffer on day of use.

Note: Allow HRP solution to warm up to room temperature before making dilutions.

Note: The Flunixin ELISA bulk kit supplies the analyst with approximately 1 mL of HRP enzyme conjugate and 200 mL of the EIA buffer. These volumes are about twice the theoretical volume needed to assay the 480 wells contained on the five plates in the bulk kit. However, the analyst is cautioned to conserve the volumes of the reagents in the kit, to allow for losses due to pipetting or mistakes. Calculate the number of wells to be used on a particular day to determine the amount of HRP conjugate to dilute. Instruction supplied with the bulk kit suggest that the analyst prepare sufficient HRP conjugate: EIA buffer solution to assay four extra wells (if using a single channel micropipettor) beyond the number of wells needed for sample extracts for the day, to account for loss by pipetting. If using a multi-channel micropipettor, allow for extra wells).

For example, if using all the wells on a 96-well plate, prepare sufficient volume of HRP conjugate: EIA buffer solution to assay 100 wells.

- e. Wash buffer, dilution – Transfer 50 mL of the concentrated wash buffer (supplied with ELISA kit at 10x concentration) into a 500-mL volumetric flask, and dilute to mark with HPLC-grade water. This solution is assigned the same expiration date as the concentrated wash buffer and shall be stored at 2-8°C.

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 5 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

**D. STANDARDS**

1. Flunixin Analytical Standard
  - a. Standard Name: Flunixin N-methyl glucamine salt (Flunixin-NMG).
  - b. Sources: Schering-Plough Corporation, P.O. Box 3182, Union, NJ. 07083-1982, Cat. No. SCH 14714, or U.S Pharmacopeia, 12601 Twinbrook Parkway, Rockville, MD 20852, Cat. No, 27460-7.
  - c. CAS Registry No: 42461-84-7 (NMG salt) and 38677-85-9 (Acid).
  - d. Chemical formula:  $C_{14}H_{11}F_3N_2O_2$  (Acid),  $C_{21}H_{28}F_3N_3O_7$  (NMG Salt).
  - e. MW: 296 Acid and 491 (NMG Salt).
  - f. Storage: Per manufacturer's instructions.
  
2. Preparation of Standard Solutions
  - a. Stock Solution A (500  $\mu\text{g}/\text{mL}$  free acid): - Weigh an equivalent of 50.0 mg of Flunixin-free acid (approximately 83.0 mg of Flunixin-NMG analytical standard; see formula below) into a 100-mL volumetric flask, dissolve the material, and dilute to the mark with methanol. This standard is stable for 6 months when stored at  $< -10^\circ\text{C}$ .  
$$\text{Weight of Flunixin NMG (mg)} = [(50.0 \times 491) / \text{purity}] / 296.$$

For a purity of 99.0%, substitute 0.990 into the above equation.
  - b. Stock Solution B (50  $\mu\text{g}/\text{mL}$  free acid): - Pipet 10 mL of stock solution A into a 100-mL volumetric flask and dilute to the mark with methanol. This standard is stable for 6 months when stored at  $< -10^\circ\text{C}$ .
  - c. Working standard solutions:
    - i. 50 ng/mL (50 ppb): - Pipet 100  $\mu\text{L}$  of stock solution B into a 100-mL volumetric flask and dilute to mark with HPLC-grade water. Mix thoroughly. This standard is stable for one month when stored at  $2-8^\circ\text{C}$ .
    - ii. 120 ng/mL - Pipet 240  $\mu\text{L}$  of stock solution B into a 100 mL volumetric flask and dilute to mark with HPLC-grade water. Mix thoroughly. This standard is stable for one month when stored at  $2-8^\circ\text{C}$ .

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 6 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

- d. External Standards. Prepare fresh daily
  - i. 15 ng/mL (Equivalent to 100 ppb in sample). Add 100  $\mu$ L 120 ng/mL working standard (D.2.c.ii) to 700  $\mu$ L extraction buffer in an appropriate vial and mix.
  - ii. 7.5 ng/mL (Equivalent to 50 ppb in sample). Mix equal volumes of 15 ng/mL (D.2.d.i) standard and extraction buffer in an appropriate vial.

**E. SAMPLE PREPARATION**

Samples of liver and muscle are processed until homogeneous. All samples are stored refrigerated (short term/overnight) or frozen (long term) until analyzed.

**F. ANALYTICAL PROCEDURE**

1. Sample Set Requirements

- a. Each sample set will contain a negative control and a positive control for the appropriate tissue type being analyzed plus any number of samples that can be accommodated in a set.

2. Control sample preparation

- a. Weigh  $1.0 \pm 0.10$  g portion of blank tissue into two 50-mL centrifuge tubes. (A blank is previously analyzed tissue found to contain no Flunixin residue and exhibit minimal background activity). Use one tube as negative control.

- a. Positive Control

- i. Liver (50 ng/g)

- For fortification, add 1-mL of 50 ng/mL working standard (D.2.c.i) to the remaining tube.

- ii. Muscle (10 ng/g)

- For fortification, add 200  $\mu$ L of 50 ng/mL working standard (D.2.c.i) to the remaining tube.

3. Extraction Procedure

- a. Weigh a  $1.0 \pm 0.10$  g portion of each frozen tissue into a 50 mL centrifuge tube. Allow to thaw.

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 7 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

i. Negative Control and Samples

Add 6 mL of phosphate extraction buffer (C.2.a) to the negative tube and sample tubes.

ii. Positive Controls

Liver: Add 5 mL of phosphate extraction buffer (C.2.a) to the positive liver tissue control tube, which already contains 1 mL of working standard.

Muscle: Add 5.8 mL of phosphate extraction buffer (C.2.a) to the positive muscle tissue control tube, which already contains 200 µL of working standard.

- b. Mix contents of tube on vortex mixer for 5 seconds.
- c. Shake tube vigorously for 5 min.
- d. Sonicate tube for 1 min.
- e. Centrifuge at 15,000 RPM at 4°C for 15 min.
- f. Keep tube refrigerated until ready for ELISA.

4. ELISA

- a. Apply 20 µL portions of control tissue and sample extracts in duplicate to the ELISA plate.

Note: An estimation of well to well repeatability is necessary to properly set a cutoff based on the response of the positive control (50 ppb for liver and 10 ppb for muscle). This may be estimated by calculating the standard deviation of the normalized absolute differences between all duplicate wells of replicate controls.

Distribute randomly the fortified control extracts in 20 µL portions into 3 pairs of wells over the area of the ELISA plate utilized.

- b. Following ELISA kit directions, add 180 µL diluted drug-HRP Conjugate Solution (HRP) to each well. (The HRP solution and EIA buffer are supplied in the test kit and are diluted 1:180 HRP: EIA or per Kit Instruction). After addition of the HRP: EIA solution to each sample, mix the solutions in plate wells by gently vibrating the plate on a flat surface.
- c. Cover plate to avoid possible dust/dirt contamination.
- d. Incubate for one hour at room temperature. Shake gently at least twice during the incubation period.
- e. Invert the plate after the incubation period, to remove matrix solutions.
- f. Wash the wells 3 times with (300 µL/time) with diluted washing buffer (diluted 1:10 with water and supplied with the test kit). Tamp the inverted plate on a paper towel between washings.

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01		Page 8 of 12
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

- g. Add 150 µL K-Blue substrate (another reagent supplied in the test kit) to each well. Allow the reaction to proceed for 30 minutes with intermittent gentle shaking of plate, especially before taking an optical density (absorbance) reading.
- h. Read results at 630nm or 650nm on a Biotec Autoreader ELx 808. The optimum absorbance value is considered to be between 1.0 to 2.0 for the negative control read against air.

**G. CALCULATIONS**

1. Evaluate sample results based on absorbance values for fortified control tissue. Flunixin concentration is inversely related to intensities of blue color.
2. Calculate the mean and standard deviation (SD) for the absorbance readings of the six positive control replicates (see F.4.a). Use these to determine a decision level (DL) using the formula  $DL = Mean + 3*SD$ . Identify sample as positive if its average absorbance is less than the decision level.

**H. HAZARD ANALYSIS**

1. Required Protective Equipment — Lab coat, safety glasses, and gloves.
2. Hazards

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Primary standard	Causes skin irritation and possible pulmonary edema after prolonged inhalation of this compound	Handle primary standard in a fume hood. Use safety glasses and gloves when handling the primary standard.
Methanol	This solvent may be flammable and may produce toxic effects to skin, eyes and the respiratory system	Use reagents in an efficient fume hood away from all electrical devices and open flames.
Acid and base	Corrosive	Wear gloves and safety glasses.

3. Disposal Procedures

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
-----------------------	---------------	------------------------------------



United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 9 of 12
Title: Screening of Flunixin Residues by ELISA	
Revision 01	Replaces: CLG-FLX3.00
Effective: 03/10/2011	

Methanol	Can cause burns on skin and eye injury	Collect waste in tightly sealed container and store away from non-compatibles in a cool, well ventilated, flammable liquid storage area/cabinet for disposal in accordance with local, state, and Federal regulations.
Acid and base	See Above	Neutralize the acid or base for disposal down the drain in accordance with local, State and Federal regulations.
Used plate(s)	None	Dispose in accordance with local, State and Federal regulations.

**I. QUALITY ASSURANCE PLAN**

1. Performance Standard

<i>Analyte</i>	<i>Analytical Range</i>	<i>Acceptable Recovery</i>	<i>Acceptable Repeatability (CV)</i>
Flunixin	≥ 50 ppb for liver ≥ 10 ppb for muscle	NA	NA

- a. Assess the following criteria:
  - i. The negative control tissue absorbance should be between 1.0 to 2.0 absorbance units vs. air.
  - ii. Absorbance must decrease with increasing Flunixin concentration for all positive controls.
  - iii. No false positives for the blank (negative control).

2. Critical Control Points and Specifications

<i>Record</i>	<i>Acceptable Control</i>
---------------	---------------------------

---

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 10 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

- a. F.3.a: Sample Weight                      Weigh  $1.00 \pm 0.10$  g.
- b. F.3.g: Extract storage condition        Store in refrigerator until ready for ELISA kit.
- c. C.2.d. and F.4.b. Horseradish        Prepare on the day of use. Dilute properly.  
peroxidase dilution

3. Readiness To Perform

a. Familiarization:

Note: Phase I and Phase II may be performed concurrently.

- i. Phase I: When possible, the analyst being trained may observe the method being performed by an analyst who has previously qualified in the method. Run duplicate external standard curves on 3 working days, which will include the following:
  - (a) 0 ppb (1 mL buffer solution).
  - (b) 50 ppb - (Working Standard (D.2.d.ii) diluted 1:7 with buffer).
  - (c) 100 ppb - (Working Standard (D.2.d.i) diluted 1:7 with buffer).
- ii. Phase II: Fortified samples- 3 sets of 10 blank and 10 fortified (50 ppb for liver or 10 ppb for muscle) samples over a period of 3 working days.
- iii. Phase III. Check samples for analyst accreditation.
  - (a) The analyst is to receive 30 samples. The sample fortifications, including the number of blanks, are to be blind to the analyst. At least 10 of the 30 samples should be blank and the rest spiked at or greater than the appropriate positive control level. These samples will be prepared and provided by Supervisor or Quality Assurance Manager (QAM). All samples should be analyzed in duplicate wells. The samples must be randomized throughout the set.
  - (b) An external curve must be generated using 0, 50 and 100 ppb, respectively to help monitor plate acceptability.
  - (c) Report analytical findings to the QAM (Quality Assurance Manager) and Supervisor.
  - (d) Authorization from QAM and Supervisor is required to commence official analysis.

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 11 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

- b. Acceptability criteria:  
Refer to section I.1
- 4. Intralaboratory Check Samples
  - a. System, minimum contents:
    - i. Frequency: Once per week per analyst when samples analyzed.
    - ii. Records are to be maintained.
  - b. Acceptability criteria:  
Refer to section I.1  
If unacceptable values are obtained, then:
    - i. Stop all official analyses by that analyst for this method.
    - ii. Take corrective action.
- 5. Sample Acceptability and Stability:
  - a. Matrix: Bovine liver and muscle.
  - b. Sample receipt size: Minimum 20 – 25 g
  - c. Condition upon Receipt: Cool or frozen
  - d. Sample storage:
    - i. Time: 2 months
    - ii. Condition: Frozen stored at < -10°C
- 6. Sample Set:
  - a. Each sample set contains
    - i. Negative Control
    - ii. Positive Control
    - iii. Samples up to 39
- 7. Minimum Level of Applicability (MLA): 50 ppb for liver and 10 ppb for muscle .

**J. WORKSHEET**  
[RESERVED]

United States Department of Agriculture  
Food Safety and Inspection Service, Office of Public Health Science

CLG-FLX3.01	Page 12 of 12	
Title: Screening of Flunixin Residues by ELISA		
Revision 01	Replaces: CLG-FLX3.00	Effective: 03/10/2011

**K. REFERENCES**

Susan B. Clark, Sherri B. Turnipseed, Gene J. Nandrea, Mark R. Madson, and Emma R. Singleton, *Identification and Confirmation of Flunixin Meglumine and Phenylbutazone Residues in Animal Kidney by ELISA Screening and Liquid Chromatography Mass Spectrometry*, US FDA Laboratory Information Bulletin, #4246, May 2001.

**L. APPROVALS AND AUTHORITIES**

1. Approvals on file.
2. Issuing Authority: Director, Laboratory Quality Assurance Division.